

Package ‘SCArray.sat’

July 2, 2025

Type Package

Title Large-scale single-cell RNA-seq data analysis using GDS files and Seurat

Version 1.9.0

Date 2025-03-24

Depends methods, SCArray (>= 1.13.1), SeuratObject (>= 5.0), Seurat (>= 5.0)

Imports S4Vectors, utils, stats, BiocGenerics, BiocParallel, gdsfmt, DelayedArray, BiocSingular, SummarizedExperiment, Matrix

Suggests future, RUnit, knitr, markdown, rmarkdown, BiocStyle

Description Extends the Seurat classes and functions to support Genomic Data Structure (GDS) files as a DelayedArray backend for data representation. It relies on the implementation of GDS-based DelayedMatrix in the SCArray package to represent single cell RNA-seq data. The common optimized algorithms leveraging GDS-based and single cell-specific DelayedMatrix (SC_GDSMatrix) are implemented in the SCArray package. SCArray.sat introduces a new SCArrayAssay class (derived from the Seurat Assay), which wraps raw counts, normalized expressions and scaled data matrix based on GDS-specific DelayedMatrix. It is designed to integrate seamlessly with the Seurat package to provide common data analysis in the SeuratObject-based workflow. Compared with Seurat, SCArray.sat significantly reduces the memory usage without downsampling and can be applied to very large datasets.

License GPL-3

VignetteBuilder knitr

BugReports <https://github.com/AbbVie-ComputationalGenomics/SCArray/issues>

biocViews DataRepresentation, DataImport, SingleCell, RNASeq

git_url <https://git.bioconductor.org/packages/SCArray.sat>

git_branch devel

git_last_commit cb8ac58

git_last_commit_date 2025-04-15

Repository Bioconductor 3.22

Date/Publication 2025-07-01

Author Xiuwen Zheng [aut, cre] (ORCID:
<https://orcid.org/0000-0002-1390-0708>),
 Seurat contributors [ctb] (for the classes and methods defined in
 Seurat)

Maintainer Xiuwen Zheng <xiuwen.zheng@abbvie.com>

Contents

SCArray.sat-package	2
CreateAssayObject2	3
NormalizeData	4
RunPCA	5
ScaleData	6
SCArrayAssay-class	8
SCArrayAssay-methods	9
scGetFiles	10
scMemory	11
scNewAssayGDS	12
scNewSeuratGDS	13
Seurat_g-class	14
Index	15

SCArray.sat-package	<i>Large-scale single-cell RNA-seq data analysis using GDS files and Seurat</i>
---------------------	---

Description

The package extends the Seurat classes and functions to support GDS files as a DelayedArray backend for data representation. It introduces a new SCArrayAssay class (derived from the Seurat Assay), which wraps raw counts, normalized expressions and scaled data matrix based on DelayedMatrix. It is designed to integrate seamlessly with the SeuratObject and Seurat packages to provide common data analysis, with the optimized algorithms for GDS data files.

Details

Package: SCArray.sat
 Type: Package
 License: GPL version 3

Author(s)

Xiuwen Zheng

CreateAssayObject2	Create an Assay object
--------------------	------------------------

Description

Create an `SCArrayAssay` (inherited from `Assay`) object from counts or prenormalized data.

Usage

```
CreateAssayObject2(counts, data, min.cells=0, min.features=0,  
  key=NULL, check.matrix=FALSE, ...)
```

Arguments

<code>counts</code>	Unnormalized raw counts (matrix, <code>dgCMatrix</code> or <code>DelayedMatrix</code>)
<code>data</code>	Prenormalized data (matrix, <code>dgCMatrix</code> or <code>DelayedMatrix</code>)
<code>min.cells</code>	if > 0, a lower cutoff for filtering cells
<code>min.features</code>	if > 0, a lower cutoff for filtering features
<code>check.matrix</code>	Check counts matrix for NA, NaN, Inf, and non-integer values
<code>key</code>	Key name for the assay
<code>...</code>	Arguments passed to as.sparse when counts or data is matrix or <code>dgCMatrix</code>

Details

Similar to `SeuratObject::CreateAssayObject()`, except allowing `DelayedMatrix` counts or data, and returning a `SCArrayAssay` object. `counts` and `data` should not be provided at the same time.

Value

Return an instance of [SCArrayAssay](#).

Author(s)

Xiuwen Zheng

See Also

[CreateAssayObject](#), [scGetFiles](#)

Examples

```
fn <- system.file("extdata", "example.gds", package="SCArray")  
  
x <- scArray(fn, "counts")  
colnames(x) <- paste0("c", 1:ncol(x))  
rownames(x) <- paste0("g", 1:nrow(x))  
x  
  
a <- CreateAssayObject2(x)  
a
```

```

scGetFiles(x)
scGetFiles(a)

remove(x, a)

```

NormalizeData

Normalize Count Data

Description

Normalizes the count data in the Seurat assay.

Usage

```

# NormalizeData(object, ...)
## S3 method for class 'SC_GDSMatrix'
NormalizeData(object,
  normalization.method="LogNormalize", scale.factor=1e4, margin=1,
  verbose=TRUE, ...)

```

Arguments

<code>object</code>	input R object (e.g., a <code>SC_GDSMatrix</code> object)
<code>normalization.method</code>	"LogNormalize", "CLR" or "RC"; see <code>NormalizeData.Seurat</code> for more details
<code>scale.factor</code>	the scale factor for cell-level normalization
<code>margin</code>	only applicable when <code>normalization.method="CLR"</code> , normalize across features (<code>margin=1</code>) or cells (<code>margin=2</code>)
<code>verbose</code>	if TRUE, show information
<code>...</code>	additional arguments passed to specific methods

Details

`NormalizeData()` does not store the normalized data in a GDS file, since the calculation is "delayed" until it is needed.

Value

Returns a `SC_GDSMatrix` matrix.

Author(s)

Xiuwen Zheng

See Also

[NormalizeData](#)

Examples

```
fn <- system.file("extdata", "example.gds", package="SCArray")

d <- scNewSeuratGDS(fn)
d
d <- NormalizeData(d)

remove(a, d)
```

RunPCA

*Run PCA***Description**

Performs PCA on a Seurat SCArrayAssay or a DelayedMatrix object.

Usage

```
# RunPCA(object, ...)
## S3 method for class 'SCArrayAssay'
RunPCA(object, assay=NULL, features=NULL, npcs=50,
        rev.pca=FALSE, weight.by.var=TRUE, verbose=TRUE, ndims.print=1:5,
        nfeatures.print=30, reduction.key="PC_", seed.use=42, ...)
## S3 method for class 'SC_GDSMatrix'
RunPCA(object, assay=NULL, npcs=50, rev.pca=FALSE,
        weight.by.var=TRUE, verbose=TRUE, ndims.print=1:5,
        nfeatures.print=30, reduction.key="PC_", seed.use=42, approx=TRUE,
        BPPARAM, ...)
```

Arguments

object	input R object (e.g., a SCArrayAssay object)
assay	NULL for using the active assay, or an assay name
features	if NULL, PCA will be run on the scaled data; otherwise, features to compute PCA on
npcs	# of top PCs to be calculated
rev.pca	By default (FALSE), perform PCA on the cell x gene matrix; otherwise, compute it on gene x cell matrix
weight.by.var	if TRUE, weight the cell embeddings (when rev.pca=FALSE) or the the gene loadings (when rev.pca=TRUE) by the variance of each PC
verbose	if TRUE, show information
ndims.print	which PCs to print genes for
nfeatures.print	# of genes to print for each PC
reduction.key	dimensional reduction key
seed.use	a random seed; or NULL for not setting a seed internally
approx	if TRUE, use IrlbaSVD; otherwise use ExactSVD
BPPARAM	NULL for non-parallel execution, or a BiocParallelParam object for parallelization; if it is missing, getAutoBPPARAM() will be used
...	additional arguments passed to specific methods

Details

`RunPCA()` computes the covariance matrix of genes (if # of genes $\leq$ # of cells) or the cell covariance matrix for the PCA calculation, which can reduce the times of accessing on-disk data.

Value

Return a data frame for reduction data (via `CreateDimReducObject`).

Author(s)

Xiuwen Zheng

See Also

[RunPCA](#), [CreateDimReducObject](#), [BiocParallelParam](#), [getAutoBPPARAM](#)

Examples

```
fn <- system.file("extdata", "example.gds", package="SCArray")

d <- scNewSeuratGDS(fn)

d <- NormalizeData(d)
d <- FindVariableFeatures(d, nfeatures=250)
d <- ScaleData(d)

d <- RunPCA(d, ndims.print=1:2)
DimPlot(d, reduction="pca")

remove(a, d)
```

ScaleData

Scale and Center the Data

Description

Scales and centers features or residuals in the dataset.

Usage

```
# ScaleData(object, ...)
## S3 method for class 'SC_GDSMatrix'
ScaleData(object, features=NULL, vars.to.regress=NULL,
  latent.data=NULL, split.by=NULL, model.use='linear', use.umi=FALSE,
  do.scale=TRUE, do.center=TRUE, scale.max=10, block.size=1000,
  min.cells.to.block=3000, verbose=TRUE, use_gds=TRUE, rm_tmpfile=TRUE, ...)
```

Arguments

<code>object</code>	input R object (e.g., a <code>SC_GDSMatrix</code> object)
<code>features</code>	if <code>NULL</code> , to use the variable features (found via <code>FindVariableFeatures()</code>); or features names to scale/center
<code>vars.to.regress</code>	<code>NULL</code> or variable names to regress out
<code>latent.data</code>	<code>NULL</code> or a <code>data.frame</code> to regress out the covariates
<code>split.by</code>	variable name in the metadata, or a vector or factor defining grouping of cells
<code>model.use</code>	regression model: "linear" (default), "poisson" or "negbinom"
<code>use.umi</code>	only applicable when the covariates are given in <code>vars.to.regress</code> for regression; default is <code>FALSE</code> for linear regression, <code>TRUE</code> for negbinom and poisson models
<code>do.scale</code>	if <code>TRUE</code> , scale the data
<code>do.center</code>	if <code>TRUE</code> , center the data
<code>scale.max</code>	max value in the resulting scaled data; see <code>ScaleData.Seurat</code> for more details
<code>block.size</code>	not used
<code>min.cells.to.block</code>	not used
<code>verbose</code>	if <code>TRUE</code> , show information
<code>use_gds</code>	if <code>TRUE</code> , use <code>SC_GDSMatrix</code> for the scaled data; if <code>FALSE</code> , to use a dense in-memory scaled matrix; or a GDS file name for storing the resulting data matrix; see details
<code>rm_tmpfile</code>	if <code>TRUE</code> , remove any temporary GDS file after the calculation; the temporary file will be created when calculating the residuals
<code>...</code>	additional arguments passed to specific methods

Details

`ScaleData()` stores the scaled data in a GDS file when `use_gds=TRUE` or an output GDS file name is given via `use_gds`. When `vars.to.regress` and `split.by` are both `NULL`, an output GDS file is not needed, since the resulting `DelayedMatrix` can be represented as common operations on the count matrix. If `use_gds=TRUE`, an output file name `"_scale_data.gds"` will be used if it does not exist, or `"_scale_data2.gds"` (if not exists), `"_scale_data3.gds"` and so on. If `use_gds` is an output file name, the resulting data matrix will be saved to a GDS file. When `vars.to.regress` are given, a temporary GDS file (e.g., `"_temp_scale_data.gds"`, `use_gds` with a prefix `"_temp"`) will be created to store the residuals before scaling. This temporary file will be deleted after the calculation when `rm_tmpfile=TRUE`.

Value

Returns a `SC_GDSMatrix` matrix if `use_gds=TRUE` or `use_gds` is an output file name, otherwise returns an in-memory matrix.

Author(s)

Xiuwen Zheng

See Also

[ScaleData](#)

Examples

```

fn <- system.file("extdata", "example.gds", package="SCArray")

d <- scNewSeuratGDS(fn)

d <- NormalizeData(d)
d <- FindVariableFeatures(d, nfeatures=50)
d <- ScaleData(d)

GetAssayData(d, slot="scale.data") # DelayedMatrix

# scale with split.by
ss <- rep(c(TRUE, FALSE), length.out=ncol(d))
d <- ScaleData(d, split.by=ss)

fn <- scGetFiles(d)
fn[2L] # the file name storing scaled data

remove(a, d)
unlink(fn[grepl("^_scale", fn)], force=TRUE)

```

SCArrayAssay-class	<i>GDS-specific Assay Class</i>
--------------------	---------------------------------

Description

The SCArrayAssay class extends the Assay class of Seurat with the new slots counts2, data2 and scale.data2 replacing counts, data and scale.data.

Slots

counts2 Unnormalized raw counts (dgCMatrx or SC_GDSMatrix), replacing Assay@counts
 data2 Normalized expression data (dgCMatrx or SC_GDSMatrix), replacing Assay@data
 scale.data2 Scaled expression data (NULL, matrix or SC_GDSMatrix), replacing Assay@scale.data

Author(s)

Xiuwen Zheng

See Also

[Assay-class](#), [Seurat_g-class](#), [GetAssayData](#), [SetAssayData](#)

Description

Gets and sets data in the Seurat Assay object.

Usage

```
## S3 method for class 'SCArrayAssay'
GetAssayData(object,
  slot=c("data", "scale.data", "counts"), ...)
## S3 method for class 'SCArrayAssay'
SetAssayData(object, layer, new.data,
  slot=c('data', 'scale.data', 'counts'), ...)

## S3 method for class 'SCArrayAssay'
subset(x, cells=NULL, features=NULL, ...)
```

Arguments

object, x	a SCArrayAssay object (inherited from SeuratObject::Assay)
layer	layer
new.data	a new data matrix (dgCMatrx or SC_GDSMatrix)
slot	data matrix in the Assay object, "data" is used by default
cells	names or indices for selected cells
features	names or indices for selected features
...	further arguments to be passed to or from other methods

Value

Return a data matrix or an instance of [SCArrayAssay](#).

Author(s)

Xiuwen Zheng

See Also

[Assay](#)

scGetFiles	<i>File names for on-disk backend</i>
------------	---------------------------------------

Description

Get a list of file names for DelayedArray with an on-disk backend.

Usage

```
scGetFiles(object, ...)  
## S4 method for signature 'Assay'  
scGetFiles(object, ...)  
## S4 method for signature 'SCArrayAssay'  
scGetFiles(object, ...)  
## S4 method for signature 'Seurat'  
scGetFiles(object, ...)
```

Arguments

object	input R object (e.g., a Seurat object)
...	additional arguments passed to specific methods

Value

Return a character vector storing file names.

Author(s)

Xiuwen Zheng

See Also

[scGetFiles](#)

Examples

```
fn <- system.file("extdata", "example.gds", package="SCArray")  
  
a <- scNewAssayGDS(fn)  
d <- Seurat::CreateSeuratObject(a)  
  
scGetFiles(a)  
scGetFiles(d)  
  
remove(a, d)
```

`scMemory`*Load Data to Memory*

Description

Loads the internal data to memory for any on-disk object.

Usage

```
scMemory(x, ...)  
## S4 method for signature 'SCArrayAssay'  
scMemory(x, slot=NULL, ...)  
## S4 method for signature 'Seurat'  
scMemory(x, assay=NULL, slot=NULL, ...)
```

Arguments

<code>x</code>	input R object (e.g., a Seurat object)
<code>assay</code>	NULL for using the active assay, or a list of assay names
<code>slot</code>	NULL for all "counts", "data" and "scale.data"; or a character vector including "counts", "data" or "scale.data"; see details
<code>...</code>	additional arguments passed to specific methods

Details

If `slot=NULL`, return a Assay object instead of SCArrayAssay object, so it can downgrade a SCArrayAssay object to a Assay object.

Value

Return an object (it maybe a different type from `class(x)`).

Author(s)

Xiuwen Zheng

See Also

[scMemory](#)

Examples

```
fn <- system.file("extdata", "example.gds", package="SCArray")  
  
d1 <- scNewSeuratGDS(fn)  
is(GetAssay(d1))  
  
d2 <- scMemory(d1)  
is(GetAssay(d2))  
  
remove(a, d1, d2)
```

scNewAssayGDS

Create Assay Object

Description

Creates a new Seurat Assay object (SCArrayAssay) from a GDS file.

Usage

```
scNewAssayGDS(gdsfile, name="counts", key="rna_", row_data=TRUE, check=TRUE,
              verbose=TRUE)
```

Arguments

gdsfile	a file name for the GDS file, or a SCArrayFileClass object
name	characters for the name of data matrix in the GDS file; if NA_character_, to use the first assay
key	characters for the Assay key
row_data	if TRUE, add rowData() to the feature-level meta data of the Seurat Assay
check	if TRUE, check the feature names
verbose	if TRUE, show information

Value

Return an instance of [SCArrayAssay](#).

Author(s)

Xiuwen Zheng

See Also

[SCArrayAssay](#), [SCArrayFileClass](#), [scExperiment](#), [scNewSeuratGDS](#)

Examples

```
# raw count data in a GDS file
fn <- system.file("extdata", "example.gds", package="SCArray")

a <- scNewAssayGDS(fn)
a
class(a)

d <- Seurat::CreateSeuratObject(a)
d

rm(a, d)
```

scNewSeuratGDS	<i>Create Seurat Object</i>
----------------	-----------------------------

Description

Creates a new Seurat object from a GDS file.

Usage

```
scNewSeuratGDS(gdsfile, assay.name=NULL, key=c(counts="rna_"), row_data=TRUE,
  col_data=TRUE, check=TRUE, verbose=TRUE)
```

Arguments

gdsfile	a file name for the GDS file, or a SCArrayFileClass object
assay.name	characters for the name of data matrix in the GDS file; if NULL, to use all of the assays
key	a character vector for an assay key map, where its names are the GDS node names
row_data	if TRUE, add rowData() to the feature-level meta data of the first Seurat Assay
col_data	if TRUE, add colData() to the cell-level meta data of the Seurat object
check	if TRUE, check the feature names
verbose	if TRUE, show information

Details

"counts" must be in the input GDS file and it is used as the raw count data in the active Seurat assay. If "logcounts" exists, it is used as normalized data associated with "counts". If there are other data matrices in the GDS file, they will be added to the assay list.

Value

Return an instance of [Seurat](#).

Author(s)

Xiuwen Zheng

See Also

[SCArrayAssay](#), [SCArrayFileClass](#), [scExperiment](#), [scNewAssayGDS](#)

Examples

```
# raw count data in a GDS file
fn <- system.file("extdata", "example.gds", package="SCArray")

d <- scNewSeuratGDS(fn)
d
class(d)

rm(d)
```

Seurat_g-class	<i>GDS-based Seurat Class</i>
----------------	-------------------------------

Description

The Seurat_g class inherits directly from "Seurat".

Slots

```
setClass("Seurat_g", contains="Seurat")
```

Author(s)

Xiuwen Zheng

See Also

[SCArrayAssay-class](#)

Index

* GDS

CreateAssayObject2, [3](#)
NormalizedData, [4](#)
RunPCA, [5](#)
ScaleData, [6](#)
SCArray.sat-package, [2](#)
SCArrayAssay-class, [8](#)
SCArrayAssay-methods, [9](#)
scGetFiles, [10](#)
scMemory, [11](#)
scNewAssayGDS, [12](#)
scNewSeuratGDS, [13](#)
Seurat_g-class, [14](#)

* SingleCell

CreateAssayObject2, [3](#)
SCArray.sat-package, [2](#)
SCArrayAssay-class, [8](#)
SCArrayAssay-methods, [9](#)
scNewAssayGDS, [12](#)
scNewSeuratGDS, [13](#)
Seurat_g-class, [14](#)

* methods

NormalizedData, [4](#)
RunPCA, [5](#)
ScaleData, [6](#)
scGetFiles, [10](#)
scMemory, [11](#)

as.sparse, [3](#)

Assay, [9](#)

BiocParallelParam, [6](#)

CreateAssayObject, [3](#)
CreateAssayObject2, [3](#)
CreateDimReducObject, [6](#)

GetAssayData, [8](#)
GetAssayData (SCArrayAssay-methods), [9](#)
getAutoBPPARAM, [6](#)

NormalizedData, [4](#), [4](#)

RunPCA, [5](#), [6](#)

ScaleData, [6](#), [7](#)

SCArray.sat (SCArray.sat-package), [2](#)
SCArray.sat-package, [2](#)
SCArrayAssay, [3](#), [9](#), [12](#), [13](#)
SCArrayAssay (SCArrayAssay-class), [8](#)
SCArrayAssay-class, [8](#)
SCArrayAssay-methods, [9](#)
SCArrayFileClass, [12](#), [13](#)
scExperiment, [12](#), [13](#)
scGetFiles, [3](#), [10](#), [10](#)
scGetFiles, Assay-method (scGetFiles), [10](#)
scGetFiles, SCArrayAssay-method
(scGetFiles), [10](#)
scGetFiles, Seurat-method (scGetFiles),
[10](#)
scMemory, [11](#), [11](#)
scMemory, SCArrayAssay-method
(scMemory), [11](#)
scMemory, Seurat-method (scMemory), [11](#)
scMemory, Seurat_g-method (scMemory), [11](#)
scNewAssayGDS, [12](#), [13](#)
scNewSeuratGDS, [12](#), [13](#)
SetAssayData, [8](#)
SetAssayData (SCArrayAssay-methods), [9](#)
Seurat, [13](#)
Seurat_g (Seurat_g-class), [14](#)
Seurat_g-class, [14](#)
subset (SCArrayAssay-methods), [9](#)